

Wearable Screen-Printed DC System for Iontophoretic Sweat Extraction

Student: Benjamin Schmid Ties

Supervisor: Prof. Luisa Petti

Collaborators: Moritz Ploner, Arvind Gurusekaran, Giulia Elli

Bachelor Thesis Defense in Electrical and Cyber-Physical Systems Engineering
22nd of July 2025

Why sweat sensing:

- Sweat as alternative to blood
- Continuous monitoring of biomarkers
- Detect diseases before symptoms
- Sweat rich in analytes, correlated to blood

Challenges:

- On-demand/continuous sweat extraction
- Sweat unavailable when
 - Stationary
 - Climate controlled environment

[1]




Commercial sweat sensing unit:
Palmsense Sensit Wearable®

Iontophoresis


- Method to extract sweat on demand
- Mild current flowing through skin
- Ion transfer through membrane
- Trigger sweating in sweat glands

State-of-the-art


- Two elliptical iontophoresis electrodes
- Large contact interface with skin
- Biosensors between electrodes
- Flexible substrate



[1] Sweat stimulation electrodes



[2] Iontophoresis electrodes



State-of-the-Art iontophoretic platforms for sweat stimulation

Research gaps and objectives

Research gaps:

- Limited studies on effect of iontophoresis current intensity on ion transport
- Unclear impact of stimulation duration on ionic permeability
- No standardized in vitro models for iontophoretic testing

Objectives:

- Develop a flexible, skin-compliant iontophoretic system
- Quantify ion transport across synthetic membranes
- Explore current intensity and duration for ion-permeability in lab-scale setting

Long-term goal:

- Enable on-demand, active sweat diagnostics via integrated wearable platforms

- Screen printed Silver/Silver Chloride (Ag/AgCl) electrodes for iontophoresis
- Thermoplastic polyurethane (TPU) substrate
- Characterization:
 - Four-point probe
 - Optical profilometer
- Strat-M membrane as substitute for skin in experiments
- Artificial sweat
- Cyclic voltammetry for monitoring ion permeation

Fabrication of electrodes

- Thermoplastic polyurethane (TPU) substrate
- Silver/Silver Chloride (Ag/AgCl) ink (90%Ag 10%AgCl)
- Screen: 80 threads/cm
- Curing in oven (100 °C, 15 minutes)





Illustration of screen-printing process




Design and printed electrodes





Fabricated electrodes

Characterization of electrodes


- Average sheet resistance $44.47 \pm 10.57 \frac{m\Omega}{sq}$
- Ag/AgCl: low sheet resistance and biocompatible
- Surface roughness, (1-10 μm at edge of print)



Four-point-probe




Optical profilometry scan




Strat-M membrane as skin alternative

- Synthetic membrane
- Model for ion permeation through skin
- 3 polymer layers modelling epidermis, dermis and hypodermis
- Used for testing skin permeation of drugs





Scanning electron microscope image of Strat-M



Corresponding layers in human skin


Experimental setup for iontophoresis

- Similar to a Franz-Diffusion cell
- Electrodes attached to Strat-M membrane
- Application of DC current
- Measure ion permeation through Strat-M
- From donor to receptor chamber




Validation of ionic permeation


- Periodic sample retrieval from receptor chamber (60 μL samples every 5 minutes)
- Analysis of ionic composition of receptor solution via cyclic voltammetry
- Cyclic voltammetry performed on a three-electrode platform (DropSense®)
 - Working electrode (WE) and counter electrode (CE): Carbon
 - Reference electrode: Ag/AgCl




Cyclic voltammetry on a three-electrode platform



Applied potential between working and reference electrode




Cyclic voltammogram



Experimental setup for ionic composition measurements


Results I

- Analyze change in area enclosed by cyclic voltammograms → Proportional to conductivity of solution
- Area increases with time ions pass through membrane



18/07/2025

Cyclic voltammograms over time




Relative change in Q_t over time for different current intensities

11

Results


- Relative change in ions transferred proportional to current
- When no current is applied, no significant ion transfer across membrane
- Maximum ion transfer at 0.3mA
- 0.3mA still far below pain threshold
- Ion transfer does not saturate over time



Conclusions and outlook

Outlook:

- Establishes foundation for future wearable, continuous, on-demand sweat sampling
- Enables development of fully integrated wearable biosensors combining sweat stimulation and analysis



Iontophoretic-equipped and microfluidics-integrated sensor for sweat analysis

Thank you for your attention Questions?

Thank you to:
Prof. Luisa Petti, Moritz Ploner, Arvind Gurusekaran, Giulia Elli

22nd of July 2025